Experiment- 4

ARKA PRAMANICK, AE21B007
Department of Aerospace Engineering
IIT Madras
Basic Aerospace Engineering lab
Bourdon Gauge Caliberation

Instructor:
Professor Dr. Manikandan Mathur
19 October,2022

1. Aim:

To verify how accurate is the pressure shown by the Bourdan's gauge.

2. Apparatus:

- Bourdon Pressure gauge:-This measures the apparent gauge pressure according to the deflection in the Bourdon's pressure tube.
- Weights:- Weights can be added in top of each other to increase the pressure in the Bourdon's gauge.
- Water:-Water is needed to carry the pressure across the plastic tube to the buordon pressure gauge.
- Holder of 0.5kg.

3. Theory:

Bourdon tube pressure gauges are used for the measurement of gauge pressures from 0.6...7,000 bar. They are classified as mechanical pressure instruments, and thus operate without any supply voltage.



Figure 1: Bourdon pressure gauge



Figure 2: Bourdon tube pressure gauge

4. Procedure:

- Read the pressure reading on the gauge before loading. This is the pressure due to the water column.
- Load the weights one by one. Take readings with 5 weights.
- We can obtain the pressure applied by each weight using the formula $P_w = \frac{weight*g}{\frac{\pi D^2}{4}}$
- Now unload the weights one by one and take pressure readings at unloading.
- Perform this experiment 5 times.




Figure 3: Bourdon pressure gauge

5. Result:

5.1 Theoretical value of pressure

$$P = \frac{mg}{\frac{\pi D^2}{4}}$$
 D=17.67mm

| S.No. | mass (in kg) | Pressure,p(in kPa) |
|-------|--------------|--------------------|
| 1 | 0.5 | 20 |
| 2 | 1.0 | 40 |
| 3 | 1.5 | 60 |
| 4 | 2.0 | 80 |
| 5 | 2.5 | 100 |

Table 1: Theoretical data for corresponding mass and pressure

5.2 Experimental value of pressure

| Mass(in kg) | P(in kPa)(loading) | P(in kPa)(unloading) | Average P(loading) | Average P(unloading) | Average P |
|-------------|---------------------|----------------------|--------------------|----------------------|-----------|
| 0.5 | 20,25,25,25,25 | 25,25,25,25,25 | 24 | 25 | 24.5 |
| 1.0 | 40,45,40,40,45 | 40,40,40,40,40 | 42 | 40 | 41 |
| 1.5 | 60,60,60,60,60 | 60,60,60,60,65 | 60 | 61 | 60.5 |
| 2.0 | 80,85,85,85,85 | 85,85,85,85,85 | 84 | 85 | 84.5 |
| 2.5 | 100,100,100,100,100 | 100,100,100,100,100 | 100 | 100 | 100 |

Table 2: Experimental data for corresponding mass and pressure

6. Graph:




Figure 4: Experimental pressure vs theoretical pressure while loading




Figure 5: Experimental pressure vs theoretical pressure while unloading




Figure 6: Average experimental pressure vs theoretical

7. Sources of error:

- Inaccurate reading due to parallax error
- Due to friction between piston and walls of tube
- Varriations in experimental conditions.
- Vibrations and overpressure.

8. Conclusion:

From experimental result I can conclude that the pressure measurement using a Bourdon Gauge is fairly accurate and advantageous for applications that allows an error of 10 percentage.

The instrument is also accurate since there is a very little variation between theoretical and experimental value.